
Balanza pesa personas
Balance pèse personnes
Weighs persons scale

RGT

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
The Health scale RGT is a scale which can both weight the body and measure the height. The designation of this product meets the OIML international recommendation R76.

The scale has the advantages of large weighing range, accuracy and high sensitivity. It is widely applicable for measuring weight and height of person in enterprises, schools, hospitals, clinics and sports departments, therefore is the basic examination device of the medical department. This scale also can be used for weighing other things.

SPECIFICATIONS


Maximum capacity	200 kg
Scale division (d)	100 g
Accuracy class	III
Range of measurement	750-2000 mm (29 1/2 - 78 3/4 in)
Scale division	5 mm (1/4 in)
Size of platform	375 x 275 mm
Overall dimension	530 x 275 x 1485 mm
Deadweight	15 kg

CAUTIONS

- This precision instrument is extremely easy to set up as all major parts are factory pre-assembled. Before starting to assemble, please unpack carefully and remove all wrappings and filters.

Carton contains:



1. Column with head
2. Platform base
3. Height measuring rod
4. Wheel module

- **NOTE:** If column bracket is bent, straighten before assembling. Slip column over column bracket on platform base. Hold column straight so that draft rod enters opening in column bracket. Screw column to column brackets (two M6 truss head screws in front of column, four in back).
- Lay down scale, with column horizontal to floor.
 1. Raise draft rod keeping out of lever's way..
 2. Push lever.
 3. When lever is in position, hook draft rod around pivot.

▪ Optional height rod installation.


1. Remove the cardboard insert from the height rod.
2. Insert the two hex screws from the hardware pack into the holes in the front of the column and tighten them with the wrench.
3. Place both height rod brackets over the two pre-installed hex head screws and pull down.
4. Use the included wrench to tighten both hex head screws. Do not overtighten the screws.

- To balance scale, move both poises to zero and turn balance screw to right or left until scale balances. For accurate weighing, place scale on level floor.

▪ Height rod operating guide

1. Before the person steps onto the scale platform, the spoon should be rotated to the horizontal position, and raised well above the person's apparent height.
2. The person may now step onto the scale platform. The spoon should be held horizontal and above the person's head.
3. Carefully lower the spoon, while keeping it horizontal, until it rests gently upon the top of the person's head. If the person is shorter than 101.5 cm (3'4"), push the latch to the right, while simultaneously pushing down on the spoon until the spoon rests horizontally upon the top of the person's head.
4. Read the height of the person as follows:
 - If the back of the spoon points to the outer height rod, then it points to the correct height.
 - If the back of the spoon points to the inner height rod, then the correct height is read at the top of the outer height rod (see "Read" arrow on the outer height rod).
5. While holding the spoon horizontally, raise the spoon above the person's head. The person may now step off of the scale platform. Hold the spoon horizontal until the person is clear of the height rod.
6. Rotate the spoon back to the vertical position and adjust the height rod back to the rest position (i.e. the spoon should be locked in place within the inner height rod and the inner rod should be at its lowest position).

Baxtran

marca propiedad de /
marque propriété de /
trade mark propriety of:

Pol. Empordà Internacional
Calle Molló, 15-16
17469 VILAMALLA
(Girona) SPAIN
T. (34) 972 527 212
F. (34) 972 527 211